

# InVite™

fruit fly lure



## Non-Toxic Gel Lure for Red & Dark Eyed Fruit Flies, Phorid Flies & Drain Flies

**InVite™ Fruit Fly Lure** gel is the quick, easy way to trap and monitor small flies at minimal cost. Packaged in a standard size bait syringe, the gel is applied directly to glueboards of any shape and size. Results will likely be best if insect monitors are used rather than rodent glue boards, since small flies can easily escape some rodent glue boards. Insect monitors trap flies best if the trap is folded so that the glue surface is facing down or the trap is placed on the underside of a bar or counter. A killing station can be created



by using InVite™ Fruit Fly Lure gel on the D-Sect glueboard and then wiping the station with LambdaStar™ 9.7% CS or FenvaStar™ EcoCap™ or non-repellant. Flies that land on the station will pick up the residual and die even if they don't get trapped.

**Label and MSDS available at [www.rockwelllabs.com](http://www.rockwelllabs.com).**

**Product Code: IFFL035**  
**one 35 g syringe with plunger**